

Cleaning and disinfecting during the COVID-19 outbreak

Clean all surfaces and equipment between every patient use.



1



Use an “S” Shape pattern to cover entire area.

2



Allow to air dry for one minute.

3



Once the wipe is dry or soiled, throw it out and replace with a new one.

4



Do not rewet the wipe or use the same wipe on multiple pieces of equipment.

DO NOT FLUSH DOWN TOILETS! DO NOT DISPOSE OF IN MACERATORS!

Please ensure that you follow the equipment manufacturer's guidelines.

GHA200266